

A new opening for 'Hibidil' users

Just a pull and two twists*
-one litre of autoclaved
ready-to-use antiseptic
in 2 seconds.



HIBIDIL
Autoclaved antiseptic solution containing
0.25% v/v Chlorhexidine Gluconate
Solution B.P. (equivalent to 0.05% w/v
Chlorhexidine Gluconate)

STERILE READY FOR USE
Single use pack for procedures requiring
large volumes of solution.
Sterile if pack is unopened.
Protect from light.

1000 ml
available now
500 ml
coming shortly

Use Antiseptic for clean skin
Administration Use without further dilution for topical application only
Warnings For external use only, not injection. Avoid contact with eyes, brain, meninges, middle ear. In aseptic procedures, disinfect container outside before opening. Discard surplus after use. Keep away from children
Side effects Idiosyncratic skin reactions have been reported but are extremely rare


Accidental ingestion Perform gastric lavage with milk, egg white, gelatine or mild soap.

PACKS	BASIC NHS PRICES
250 x 25 ml sachets	£19.95
60 x 100 ml sachets	£11.45
10 x 500 ml bottles	£11.30
10 x 1000 ml bottles	£12.30


Contracts standing orders enable us to substantially reduce these prices!

Product Licence No. PL0029/0143

Hibidil is a trade mark. 523
*Further information is available on request


Imperial Chemical Industries PLC
Pharmaceuticals Division
Alderley House, Alderley Park
Macclesfield, Cheshire





Ahead of its time

New API Z reinforces API's reputation for recognising problems and producing fast, simple, cost effective solutions. In this case, screening faecal isolates for Salmonella and Shigella which was, up until now, time consuming — but no longer.

API Z uses new and unique chromogenic enzyme substrates to provide rapid and reliable results in just two hours. Simple to set-up and reducing the need for further tests. API Z sets a new high standard for faecal screens.

api

Laboratory Products Ltd

Freepost Grafton Way
Basingstoke Hants RG22 6HY
Tel. No. Basingstoke (STD 0256) 51881